

ISGAS

Custom Calibration Standards

Our standards will give you peak performance!!!

GC Standards


- ◆ Oxygenates—ASTM D5599 & D4815
- ◆ MTBE—ASTM D5441
- ◆ Ethanol—ASTM D5501
- ◆ PONA & PIONA standards
- ◆ SIMDIS—ASTM D3710 & D2887
- ◆ Benzene Purity—ASTM D5713 & D4492
- ◆ Benzene Calibration—ASTM D3606
- ◆ Aromatics in Gasoline—ASTM D5769
- ◆ Olefins in Gasoline—ASTM D6550
- ◆ Total Aromatics in Gasoline—ASTM D5580
- ◆ Aromatics in Diesel—ASTM D5186
- ◆ Boiling Range—ASTM D3710
- ◆ Residue in LP Gas—ASTM D2158
- ◆ Hydrocarbons in Gasoline—ASTM D2427
- ◆ Aromatic Hydrocarbon—ASTM D5443
- ◆ Petroleum Naphthas—ASTM D5134
- ◆ Impurities in Ethylbenzene—UOP 755
- ◆ Boiling Range—UOP 621
- ◆ Octanes in Gasoline—UOP 690
- ◆ Hexanes in Gasoline—UOP 551
- ◆ Trace Saturates—UOP 868
- ◆ Aromatics in Hydrocarbons—UOP 744
- ◆ Hydrocarbons in Sulfolane—UOP 831
- ◆ P-Diethylbenzene in C8's—UOP 798
- ◆ Impurities in Xylenes—UOP 931
- ◆ Impurities in p-Xylene—UOP 720
- ◆ Liquid Hydrocarbons—ASTM D2597
- ◆ Natural Gas—GPA 2286
- ◆ Natural Gas—GPA 2261
- ◆ Natural Gas—GPA 2261
- ◆ Natural Gas—GPA 2177
- ◆ Natural Gas—GPA 2165
- ◆ Natural Gas—ASTM D1945

Sulfur & Nitrogen Standards

We carry a full Range of ULSD Standards


- ◆ X-Ray Standards—ASTM D2622, D4294 & D7039
- ◆ Sulfur by Radiometric Change—ASTM D4045 & D6212
- ◆ Sulfur by Electrochemical Detection—ASTM D6366 & D6428
- ◆ Sulfur by UV Fluorescence—ASTM D5453
- ◆ Sulfur by High Heat Combustion—ASTM D1552
- ◆ Sulfur by Microcoulometry—ASTM D3120, D3246 & D3961
- ◆ Sulfur in Gasoline—ASTM D6334 & D6445
- ◆ Sulfur in Natural Gas—ASTM D5504
- ◆ Sulfur in LPG/Gas—ASTM D 6667
- ◆ Sulfur by SCD—GPA 2199
- ◆ Sulfur in Petroleum Liquids—ASTM D5623
- ◆ Trace Thiophene in Benzene—ASTM D7011
- ◆ Mercaptan Sulfur—ASTM D 5723 & UOP 163
- ◆ Nitrogen by Chemiluminescence—ASTM D4629 & D5762
- ◆ Chlorides in Hydrocarbons—ASTM D4929 & D5808

- Gas Blends
- LPG/LNG Blends
- Piston Blends
- Liquid Blends
- ASTM Blends
- Custom Ampouling
- Custom Product Purification
- Constant Pressure Cylinders
 - Superior Service
 - Superior Delivery Times



ISGAS

5807 Northdale St.
Phone: 713-645-5886
Fax: 713-645-5886
E-mail: sales@isgas.us
www.isgas.us



Turn-Around Time 3-5 Days for Hydrocarbons & 5-7 Days for Sulfurs

Our chromatograms separates us from the competition